
Overview of software system blueprints

What are the functions required to display images?

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ImageAnalysis software field 1

The screenshot displays the ImageAnalysis software interface within a web browser. The browser's address bar shows the file path: `C:/Users/daisu/OneDrive/ias/index.html`. The software interface includes a left sidebar with various toolbars (File/Edit, Open, Save/Load, Copy/Paste, Comment, ChangeDimension) and a central workspace. A 'VesselSelect' dialog box is open in the center, showing tabs for Slide, Dish, WellPlate, Wafer, and Custom. The WellPlate tab is active, displaying a grid of wells with counts 6, 12, 24, 48, and 96. A red box highlights the 'Vesselの自動選択' (Automatic vessel selection) area. On the right side, a 'ViewControl' panel contains several settings: a 96wellPlate grid, a magnification dropdown (4X, 10X, 20X, 40X, 100X), a ViewMethod section with color-coded buttons (S, B, G, R, C, Y, M), a Light section with three sliders (all set to 50), a Z-Position slider, and a Timeline slider. Red annotations with arrows point to specific features: '対物レンズの紐づけ・知物レンズ情報の表示' (Objective lens assignment / Objective lens information display) points to the magnification dropdown; 'チャンネルの紐づけ、ON・OFFで表示表示の切り替え' (Channel assignment, ON/OFF display toggle) points to the ViewMethod buttons; 'Zポジションの紐づけ、表示' (Z-position assignment, display) points to the Z-Position slider; and 'タイムラインの紐づけ・表示' (Timeline assignment, display) points to the Timeline slider. A red text box at the top right explains the grid selection: 'カラム・ロウ (A1やB3などの紐づけ) 紐づいた部分は青で表示してクリックするとその画像が真ん中のカラムに表示' (Column/Row (A1 or B3 etc. assignment) The assigned area is displayed in blue, and clicking it displays the image in the center column).

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

2D Filter

【emphasis】

Low pass
High pass
Sharpening
Median
Gauss
High gauss
Local equalization
Flattening
Rank
Stain removal
Sigma
Sigma(Median)

【edge】

Sobel
Roberts
Laplacian
Variance
Sobel phase
Horizontal edge
Canny
Sobel
Gabor
Differential interference restoration

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