
Overview of software system blueprints

What are the functions required to display images?

Linking image metadata

Bioformats

Using the function of Bioformats, each file (from meta information and file name) is read, Vessel (container selection, column, row binding), objective lens, channel, Z stack, time run is read and optimized.

Enable to display images in cooperation with ui.

The flat image also reads tile image data like Map and arranges them in the correct order so that the connection can be displayed.

Views of tile data

Bioformats

leaflet

OpenCV.js

Make tile images connectable and stitchable for display.

Tile images are often 100,000x100,000, so it is necessary to reduce the weight.

ex) Automatic creation of MAP hierarchy

Views of 3D and Timelapse

OpenCV.js

OpenGL

3D image and timeline display
Channel timeline, 3D function can be linked and displayed

Image insertion and display specifications

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ImageAnalysis software field 1

The screenshot shows the ImageAnalysis software interface. On the left is a menu with options like File/Edit, Open, Save/Load, Copy/Pesto, Comment, and ChangeDimension. The main area displays a 'VesselSelect' dialog box with tabs for Slide, Dish, WellPlate, Wafer, and Custom. Below these tabs are icons for different vessel types with counts: 6, 12, 24, 48, and 96. A red box highlights the WellPlate section, with the annotation 'Vesselの自動選択' (Automatic vessel selection). On the right is a 'ViewControl' panel with a 96wellPlate grid, magnification options (4X, 10X, 20X, 40X, 100X), a ViewMethod section with color-coded buttons (S, B, G, R, C, Y, M), and sliders for Light, Z-Position, and TineLine. Red annotations point to these sections: '対物レンズの紐づけ・知物レンズ情報の表示' (Objective lens assignment, objective lens information display) points to the ViewMethod section; 'チャンネルの紐づけ、ON・OFFで表示表示の切り替え' (Channel assignment, ON/OFF display switching) points to the Light section; 'Zポジションの紐づけ、表示' (Z-position assignment, display) points to the Z-Position slider; and 'タイムラインの紐づけ・表示' (Timeline assignment, display) points to the TineLine slider. A red box highlights the grid in the ViewControl panel, with the annotation 'カラム・ロウ (A1やB3などの紐づけ) 紐づいた部分は青で表示してクリックする際の画像が真ん中のカラムに表示' (Column/Row (A1 or B3 etc. assignment) assigned parts are displayed in blue, and the image shown when clicking is the image in the center column).

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